



Filtek™

Supreme Flowable Restorative



citydentalsupplies

Now our most popular flowable is easier to inject.


3M™ Filtek™ Supreme Ultra Flowable Restorative is available in a new ergonomic syringe. Designed for improved comfort and easier injection, this new product is called 3M™ Filtek™ Supreme Flowable Restorative.

Reduce bubbles and waste

Virtually no bubbles or material “run-on” during dispensing gives you better control.

No more guessing

The blue syringe barrel indicates the remaining material volume. When your syringe is empty, only the white plunger shows.



Bendable cannula


Easy-to-bend without kinking for better access.

Patented TRUE nanotechnology

Makes the nanoparticles more uniform in size, resulting in excellent wear resistance and polished retention.

Easy to hold and inject

Triangle finger plate and plunger make it quick and easy to personalise your grip.



Key benefits

- Virtually no bubbles or material “run-on”
- Easy-to-bend cannula for improved access
- Ergonomic design is easy to hold and inject
- Versatile for a wide range of indications and techniques
- Excellent adaptation and low shrinkage stress
- Can be warmed for injection molding. Warm up to 70C (158F)³ for up to one hour, up to 25 times

A shade-matched palette

3M™ Filtek™ Supreme Flowable Restorative shades are formulated to be an identical match to 3M™ Filtek™ Universal Restorative and 3M™ Filtek™ Supreme XTE Universal Restorative. (A1, A2, A3, A3.5, A4, B1, B2, C2, D2)



A2
3M™ Filtek™
Supreme Flowable
Restorative



A2
3M™ Filtek™
Universal
Restorative



A2B
3M™ Filtek™
Supreme XTE
Universal Restorative

Clinical cases

3M™ Filtek™ Supreme Flowable Restorative offers a smooth transition from material to tooth for esthetic restorations.



Before

Photos courtesy of Dr. Gunnar Reich, Munich, Germany



Flowable placement




After

3M™ Filtek™ Supreme Flowable Restorative adapts easily and seals the margin.



Before

Photos courtesy of Dr. Giuseppe Chiodera, Brescia, Italy



Flowable placement



After

Visit 3M.com.au/FiltekFlowables



3M Oral Care

3M Australia Pty Limited
Building A, 1 Rivett Road
North Ryde, NSW 2113
1300 363 454
3M.com.au



3M New Zealand Limited
94 Apollo Drive
Rosedale, Auckland 0632
0800 808 182
3M.co.nz



Indications

(Incremental up to 2 mm)



Class III



Class V



Minimally Invasive
Dentistry (MID)¹



Base



Liner



Repair of Indirect
(small defects)



Repair²



Undercut
Block-out



Pit & Fissure
Sealant

Refer to Instructions for Use
for full list of indications.

Ordering information



3M™ Filtek™ Supreme Flowable Restorative Refills

Syringe refills include: 2 – 2 g syringes;
20 dispensing tips.

SHADE	ITEM NO.
A1	6032A1
A2	6032A2
A3	6032A3
A3.5	6032A3.5
A4	6032A4
B1	6032B1
B2	6032B2
C2	6032C2
D2	6032D2
White	6032W
Extra White	6032XW
OA3	6032AO3

ITEM NO.	PRODUCT INFO.
3700T	Filtek™ Supreme Flowable 20 Dispensing Tips Grey, 20-gauge
3700T-100	Filtek™ Supreme Flowable 100 Dispensing Tips Grey, 20-gauge

1. Including small, non stress-bearing occlusal restorations. 2. Repair of acrylic temporary and composite resin. 3. Data on file.
Please always refer to the 3M Product Instructions For Use. 3M, Elipar, Filtek and Scotchbond are trademarks of 3M. Please recycle. Printed in Australia © 3M 2021. All rights reserved. AV011477401.